

| Aim  | Parameters              | Sample                                | Protocol  | Output  | Action   | Comments  |
|--|-------------------------|---------------------------------------|---|---|--|---|
| Objectives                                 | clean                   | -                                     | clean all objectives carefully and examine  |   |  |   |
| system                                     | microscope warm-up time | beads                                 | multi-point time lapse  |   | shifts in position with time                                     |   |
| lasers                                     | laser warm up time      | power meter                           | check laser stability upon engaging   |   |  |   |
|  | laser stability         | power meter/mirror slide              | Obj: 10X. Adjust AOTF to mid gray values. <b>Detector on counting mode.</b> Maximal zoom. run time lapse for 2h, 20sec intervals. Sequential scan- 1 laser at a time  | plot Average Int±SD with time   |  | acceptable range<10% for long times. 3% for 5 min   |
| measure alignment of excitation light path | laser intensity         | power meter                           | Obj: 10X. Keep all parameters fixed (low gain, avoid saturation) and record each laser line at varied intensities. <b>Pinhole open</b>  | overlay to extract intensity fluctuations   | compare intensity fluctuations to signal fluctuations in samples | focus on a slide first then apply the same for meter  |
|  | pinhole                 | power meter/ tetraspec bead           | change pinhole from 1AU to 2AU and measure intensity (<x3)  |   | adjust pinhole position  |   |
|  | detector                | mirror slide                          | with fixed laser intensity- compare the gain required for saturation  |   |  |   |
| lasers                                     | axial resolution        | mirror slide                          | 20nm wide detection windows centered around each laser line. Focus to the slide plan, change to XZ line scan. <b>Pinhole at 1AU and 0.5AU.</b>  | line profile  |  | preparation: glue the mirror to a glass slide and cover with #1.5, elvanol  |
|  | spectral separation     | mirror slide                          | Obj: 10X ; lasers: 488, 561, 633. focus on surface, open pinhole and increase zoom to 8X. <b>Adjust AOTF to avoid saturation at all scans.</b> sampling: 256*256 . Run XyZ 500-670, 5nm intervals   | plot average image intensity as a function of wavelength  |  |   |
|  | Z galvo Stability       | mirror slide                          | Excitation:488nm. Obj:40X-63X. Pinhole 1AU and 0.5AU. Focus on slide, switch to XZT mode. Acquire Z stack at <b>30nm intervals</b> , 256*256. 2 sec interval for 10min-1h.  | plot peak position and FWHM over time   |  |   |
| objectives, laser alignment, scanner       | field illumination      | chroma slide (pink)                   | Obj: All. Main laser lines. Identify slide surface and penetrate : <b>10X:75µm ; 20X:50µm ; 40X:40µm ; 63X:30µm.</b> minimal zoom, pinhole=1AU detector on counting mode  | apply smooth images, contrast enhancement and check intensity distribution. Divide by a homogen   |  | preparation: cover with #1.5, elvanol   |
| Scanners                                   | Image Distorsion        | Grid                                  | image a square at high zoom with good sampling  | check for any deviation from square shape   | create an edge map   |   |
| Co-localization                            | chromatic shift         | 1 µm tetraspec eads                   | pinhole 1AU , maximal zoom (apply over sampling), low laser power and low gain. Acquire Z stack (10X: 1µm ; 63X: 0.15µm) at sequential scanning for several laser lines   | check shift between channels (X-Y-Z) and apply a correction factor for each direction. Project this correction to next images acquired and check shifts |  |   |
| Stage                                      | stage stability         | 1µm tetraspec beads                   | 20X air objectives, exc 633 (avoid bleaching). Pinhole 1AU, XY pix=100nm. <b>1. multi point large distances; 2. place 1 bead at different locations</b> (known distances) along FOV. Time lapse for 30 min, 20sec intervals.  | filter to reduce noise, binary mask. Compare centers with time  |  |   |
| objective chromatic abb                    | PSF                     | 100nm beads directly on #1.5, elvanol | 488 laser line (intensity to give 75%image intensity from max bit range, use LUT), low noise (reduce gain). acquire Z stack (2X size of PSF)of an isolated bead closest to the coverslip plane with different objectives. Pixel size: (NA=0.75/138*138*716)(na=0.95/87*87*208)(NA=1.1/75*75*246)(na=1.4/59*59*152). <b>Pinhole:5AU, 1AU and 0.5AU</b> | generate PSF from image.  | Check homogeneity. Make sure the entire F                        | *check if correction collar is adjusted properly. Allow the slide to sit for 30min before acq. Take several PSFs for each objective |
| transmitted alignment                      | condensor-objective     | tissue section                        | 10x, kohler alignment. Center   | even  |  | more critical for wide-field  |